

# MagCore® Nucleic Acid Extraction Kits



# MagCore® Automated Nucleic Acid Kits Selection Guide



		CatNo. 24 preps	CatNo. 36 preps	CatNo. 72 preps	CatNo. 96 preps	Whole Blood	Buffy Coat	Plasma/Serum	Cell-free body fluids	Urine	Free circulating DNA	Amniotic Fluid	Cultured Cells	Animal Tissue	Plant Tissue	Bacteria/Sputum	Swab	FFPE	Stool	Forensic Specimens	PCR Product		
NGS	<b>701</b> MagCore® NGS Auto Size-Select purification Kit		MNC-01		MNC-02																■		
	<b>101</b> MagCore® Genomic DNA Whole Blood Kit (Speedy installation)		MGB400-01		MGB400-02	■	■																
	<b>102</b> MagCore® Genomic DNA Whole Blood Kit		MGB400-03		MGB400-04	■	■																
	<b>104</b> MagCore® Genomic DNA Large Volume Whole Blood Kit (1.2 ml)				MGB1200	■																	
	<b>105</b> MagCore® Plasma DNA Extraction Kit (1.2 ml)				MPD1200			■	■		■												
	<b>106</b> MagCore® Genomic DNA Whole Blood Kit (For Genotyping)		MGB400-07		MGB400-08	■																	
Genomic DNA	<b>110</b> MagCore® Cultured Cells DNA Kit		MCC-01		MCC-02							■	■										
	<b>115</b> MagCore® Circulating DNA large volume kit (3-4 ml)	MPD4000-01			MPD4000-03			■	■		■												
	<b>301</b> MagCore® Genomic DNA Plant Kit		MGP-01		MGP-02										■								
	<b>401</b> MagCore® Genomic DNA Tissue Kit		MGT-01		MGT-02									■			■	■	■	■	■	■	
	<b>405</b> MagCore® Genomic DNA FFPE One-Step Kit		MGF-01	MGF-03															■				
	<b>406</b> MagCore® Forensic DNA Direct Kit			MFC-03																		■	
	<b>502</b> MagCore® Genomic DNA Bacterial Kit		MBB-01		MBB-02												■						
	<b>504</b> MagCore® Gut Microbiome DNA Kit		MGM-01	MGM-02																		■	
	Viral Nucleic Acids	<b>202</b> MagCore® Viral Nucleic Acid Extraction Kit (Low PCR Inhibition)		MVN400-03		MVN400-04			■	■	■												
		<b>203</b> MagCore® Viral nucleic acid extraction kit (high sensitivity)		MVN400-05		MVN400-06			■	■	■												
<b>210</b> MagCore® Viral Nucleic Acid Large Volume Extraction Kit (2.4 ml)					MVN2400			■	■	■													
<b>211</b> MagCore® Viral Nucleic Acid Large Volume Extraction Kit (1.2 ml)					MVN1200			■	■	■													
Total RNA	<b>605</b> MagCore® Total RNA FFPR One-step Kit	MRF-01		MRF-03																		■	
	<b>620</b> MagCore® miRNA Extraction Kit	MMR-01		MMR-03				■	■														
	<b>631</b> MagCore® triXact RNA Kit	MRC-01		MRX-03		■							■	■	■								